

## EMERGENCY OXYGEN KIT - INVENTORY, MAINTENANCE AND ORDER FORM

O2 KIT #	DIVING UNIT NAME	DIVING UNIT LOCATION	UNIT DIVING SUPERVISOR NAME
----------	------------------	----------------------	-----------------------------

ITEMS ISSUED	SERIAL NUMBER	ITEMS ISSUED	QTY	ITEMS ISSUED	QTY
FIRST STAGE REGULATOR	SERIAL NUMBER	REUSABLE O2 MASK		NON-REBREATHER MASK (single use)	
SECOND STAGE DEMAND VALVE REG with LP HOSE	SERIAL NUMBER / MANUFACTURER	BVM RESUSCITATOR BAG		SPARE WASHER	
OXYGEN CYLINDER #1 **	SIZE	Steel Alum	HYDROSTATIC TEST DATE	SERIAL NUMBER	
OXYGEN CYLINDER #2 **	SIZE	Steel Alum	HYDROSTATIC TEST DATE	SERIAL NUMBER	

\*\*Cylinders must meet hydrostatic testing requirements

INSPECTION / MAINTENANCE		YEAR	JAN	FEB	MAR	APR	MAY	JUN	JUL	AUG	SEP	OCT	NOV	DEC
1.	Check condition of reusable O2 mask and the BVM resuscitator bag and mask. Clean or replace if necessary.													
2.	Check condition of single use non-rebreather mask. Replace if necessary.													
3.	Verify washer is seated properly. Tighten as needed.													
4.	Check cylinder condition. Valve should be clean and dry (NO OIL)													
5.	Check hydrostatic test date. Aluminum cylinders within 5 yrs; Steel cylinders within 5 yrs or 10 yrs w/ *													
6.	Verify cylinder wrench is attached to yoke and accessible to both cylinders (when applicable).													
7.	Open cylinder valve. Check pressure gauge. Verify cylinder is full (2015 psi +/- 10%).													
8.	Briefly push button on demand valve. Valve should have high flow, no whistle, and fully stop when released.													
9.	Close cylinder valve. Check for leaks. Pressure gauge needle should not move.													
10.	Push button on demand valve to bleed off contents.													
11.	Check for corrosion. Rinse and clean as needed.													
12.	Check general condition of watertight storage case.													
		<b>INSPECTOR'S INITIALS:</b>												
<b>DEMAND VALVE TESTS:</b> PASSING PRESSURE = 32-48 mmHg	PREVIOUS 2 YEAR SERVICE DATE	JAN/FEB/MAR			APR/MAY/JUN			JUL/AUG/SEPT			OCT/NOV/DEC			
	DATE	DATE	DATE	DATE	DATE	DATE	DATE	DATE	DATE	DATE	DATE	DATE	DATE	DATE
		mmHg		mmHg		mmHg		mmHg		mmHg		mmHg		mmHg

- NOTES:
- Emergency oxygen kits shall be inspected before diving operations and monthly inspections are recommended.
  - Follow the manufactures recommendations for the demand valve resuscitator or it shall be tested every three (3) months and serviced by an authorized service center every two (2) years.

For order information, questions or service needs for this kit, contact the NOAA Dive Center at 206-526-6623.	NAME of INSPECTOR	SIGNATURE of INSPECTOR	DATE
---	-------------------	------------------------	------